

## Typical Life and Maintenance Cycles for Cutting Components of the CG03 Cutter Grinder

Cutting Component	Life/Maintenance Cycle
Precut Blade	<p>1<sup>st</sup> Usage Turn over at 4K samples                      2<sup>nd</sup> Usage Sharpen at 8K samples                      3<sup>rd</sup> Usage Turn over after 12K samples                      4<sup>th</sup> Usage Replace after 16K samples</p> <p>If sharpened properly, and the blade is not bowed in any way, then it MAY be possible to resharpen at 16K samples to give a life of 24K samples. Minimum sharpening limits to be observed as per the User's Manual.</p>
Main Blades	<p>1<sup>st</sup> Usage Sharpen at 4K samples                      2<sup>nd</sup> Usage Sharpen at 8K samples                      3<sup>rd</sup> Usage Replace after 12K samples</p> <p>If sharpened properly and if the blades are in good condition, then it MAY be possible to resharpen at 12K samples to give a life of 16K samples. Blade Retaining Bolts MUST be replaced every 4K samples. Minimum sharpening limits to be observed as per the User's Manual.</p>
Stators	<p>1<sup>st</sup> Usage Sharpen at 8K samples                      2<sup>nd</sup> Usage Replace after 16K samples</p> <p>If sharpened properly and if the stators are in good condition, then it MAY be possible to resharpen at 16K samples to give a life of 24K samples. Minimum sharpening limits to be observed as per the User's Manual.</p>
Screen Plate	<p>1<sup>st</sup> Usage Replace after 12K samples</p> <p>Because of the surface hardening treatment, the screen plate CANNOT be sharpened, reversed or surface ground.</p>
Ejector	<p>1<sup>st</sup> Usage Replace after 32K+ samples</p> <p>Because of material limitations, the ejector CANNOT be reversed or surface ground. Its life will depend on the abrasion in the system and must be judged by each mill. The top surfaces of the "wings" of the ejector may be "hardfaced" using welding techniques and remachined by competent tradespeople for extended life.</p>

**WARNING:** These are recommendations only (and **NOT guaranteed performance figures**) for clean cane under normal laboratory or mill conditions. If cane contains dirt or other abrasive materials, the life of the cutting components will be reduced. Average sample size is 3-5Kg.